



DOSIMETRY

Gantry Mount

Accessory for the myQA[®] Daily or MatriXX Resolution[™] Detectors

The Gantry Mount enables easy and stable attachment of your 2D ionization chamber array to the Linac head, enabling advanced Machine and Patient QA under any Linac angle.

The following detectors are compatible with the new Gantry Mount:

- myQA[®] Daily for your morning Machine QA
- MatriXX Resolution[™] for your advanced Machine QA as well as for Patient QA

Benefits

- Fast, accurate and easy to use locking mechanism
- Perpendicular Machine and Patient QA measurements for any gantry rotation angle.
- Comfortable fixations hold the 2D array precisely in its position
- Compatible with Varian[®] and Elekta[®] Linacs and for different SSD needs
- Upgrades of previous Gantry Mounts for StarTrack and MatriXX Evolution are available



Fig.: Gantry Mount with detector array with gantry at 0° position (left) and at 270° position (right)

Material	<ul style="list-style-type: none">• Aluminum
Features	<ul style="list-style-type: none">• Wireless and wired use of the specified array• High precision knobs for accurate positioning• Lightweight aluminum design• Compatible with RW3 build up plates

For further information or questions please [contact](#) your local IBA Dosimetry sales representative.